

# Hybrid ZnO-AuDAPT Nanoparticle Coating for Enhanced Antimicrobial Performance of Clear Orthodontic Aligners


Jelena Cheng, Mario Tsai, Nick Ojeda, Sohum Mehta, Trevor Swan

Group 2

# Introduction

- Clear aligners are vulnerable to bacterial colonization and biofilm formation, potentially causing dental complications like enamel demineralization and gingival inflammation
- Our proposed solution employs **a hybrid antimicrobial coating**, strategically combining **ZnO nanoparticles with low concentrations of AuDAPT** embedded within **a UV-crosslinked acrylic matrix**
- Aim to leverage ZnO's robust antibacterial properties alongside AuDAPT's powerful anti-biofilm capabilities

**Figure 1.** Preparation of AgNPs-PDA-ZnO-PDA-SS/AG (MPAZ) film with enhanced mechanical performance and antibacterial activity.



# Competition/Need

## Current Coatings – Limitations:

- **ZnO**: Clear but loses effectiveness after 24 hours.
- **MgO**: Antibacterial but causes visible discoloration.
- **Quaternary Ammonium**: Easy to apply but lacks long-term data.
- **AuDAPT**: Strong anti-biofilm but turns aligners purple and may leach thiols.

## Why a New Solution?

- No current option balances **clarity**, **durability**, and **safety**.
- Patients need a **transparent, long-lasting, and biocompatible** antimicrobial coating.

# Materials:

## **Zinc Oxide (ZnO) Nanoparticles**

- Properties:
  - Well-established antibacterial, antifungal, and antiviral properties
  - Inherently transparent, maintaining aesthetic quality
  - Biocompatible

## **Low-concentration AuDAPT**

- Properties:
  - Anti-biofilm and antibacterial properties, effective against drug-resistant pathogens
  - Potential to cause discoloration of the aligner
  - Possible leaching of thiol functional groups

## **UV Crosslinked Acrylic/Carbon Matrix**

- Properties
  - Anchors nanoparticles securely to the aligner surface
  - Physically entraps AuDAPT, reducing risk of leaching

# Coating Method Overview

## Low-concentration AuDAPT

- Minimizes discoloration
- Maximizes anti-biofilm & antibacterial effects

## Dip-Coating Method


- Ensures **thin, uniform** coating

## UV-Crosslinked Acrylic/Carbon Matrix


- **Anchors AuDAPT** to aligner surface
- **Enhances wear resistance** and durability
- **Reduces leaching** of gold nanoparticles & thiol groups

## Outcome:

- Clear, **long-lasting**, and **clinically safe** antimicrobial coating for aligners
  - Evaluates coating robustness and substrate bonding.



Dip coating method  
(Perovskite Photovoltaics Basic to Advanced Concepts and Implementation  
Book • 2018)



Expected process  
(UV cross-linking of unmodified DNA on glass surfaces  
10.1007/s00216-009-3045-9  
Analytical and bioanalytical chemistry)

# Coating Characterization Overview

## Optical Properties: Transparency & Band Gap

### UV-Vis Spectroscopy & Photoluminescence (PL):

- Assess transparency, color (e.g., purple hue), and band gap.
- PL identifies defect states or multiple band gaps.

## Thickness & Surface Morphology

**Alpha-step:** Measures ~100 nm thickness, uniformity.

**AFM:** Nanoscale surface details, localized thickness.

**SEM/TEM (Cross-section):** Structural validation post initial thickness scans.

## Material Uniformity

### EDX Mapping (SEM/TEM):

- Elemental distribution of ZnO.
- Detects agglomeration or uneven deposition.

## Adhesion Testing

### Tape & Scratch Tests:

- Evaluates coating robustness and substrate bonding.

# Functional Testing & Application Performance

## Crosslinking Efficiency

**FTIR:** Tracks functional group changes post UV curing.

**Raman & XPS:** Analyzes ZnO-polymer interactions, carbon bonding states.

## Antimicrobial Testing

### Confocal Microscopy + Fluorescent Labeling:

- Compares hybrid-coated, ZnO-only, AuDAPT-only, and uncoated samples.

**Ilastik Image Analysis:** Quantifies bacterial colonization.

## Biocompatibility

### Cytotoxicity Assays (MTT, Live/Dead):

- Tests on oral cells for viability, morphology, and proliferation.

## Durability & Longevity

**Mechanical:** Simulated 7-day brushing to test wear resistance.

**Chemical:** 7-day immersion in artificial saliva for leaching & stability.

# Overall Impact & Conclusion

## On Consumers...

### Health and Comfort:

- Reduced risk of oral issues like gingivitis
- Maintained transparency and comfort that consumers desire
- Reduced frequency of dental visits

## On the Environment...

### Long-lasting with minimal waste:

- UV-crosslinked coating minimizes waste

**Pollutants:** Low concentration of thiols.

## On Manufacturing...

### Dip-coating:

- Can be applied directly after the current manufacturing process without intense capital reinvestment

## Potential Hurdles and Areas for Future Research

**Regulatory Frameworks:** Nanoparticle-coated dental appliances may require robust long-term data

**Adhesion:** Acrylics and TPU may pose an issue

### Next Steps:

- Optimizing nanoparticle size if needed
- Long-term clinical trials
- Keep a keen eye on the thiol-modified AuDAPT component





# References

- [1] Gharibnavaz, M., Arash, V., Pournajaf, A., Najafi, F., Rahmati Kamel, M., & Seyedmajidi, S. (2025). Study on the antibacterial properties and optical characteristics of clear orthodontic aligners coated with zinc oxide and magnesium oxide nanoparticles. *Orthodontics & Craniofacial Research*. <https://doi.org/10.1111/ocr.12899>
- [2] Jiang, W. (2025). Effect of bioceramic fillers on antibiofilm and remineralization properties of flowable composite used as clear aligner attachments (T). University of British Columbia. Retrieved from <https://open.library.ubc.ca/collections/ubctheses/24/items/1.0448057>
- [3] Giannini, L., Galbiati, G., Tartaglia, F. C., Grecolini, M. E., Maspero, C., & Biagi, R. (2025). Orthodontic Treatment with Fixed Appliances Versus Aligners: An Experimental Study of Periodontal Aspects. *Dentistry Journal*, 13(2), 70. <https://doi.org/10.3390/dj13020070>
- [4] Thiruba Suhasini, P., Chakravarthi, S. N. C., Davis, D., Bharathi, E., & Rasiga, G. (2025). Aligner materials and their properties. *Cuestiones de Fisioterapia*, 54(4). <https://doi.org/10.48047/1yhy8771>
- [5] Gramuglia, V., Makowka, S., Tanberg, W., Zhou, A., Sharma, A., & Al-Jewair, T. (2025). Effect of quaternary ammonium-based antimicrobial coating on the mechanical properties and bacterial adhesion to Gummetal archwire. *International Journal of Dentistry*, 2025, Article 8157347, 8 pages. <https://doi.org/10.1155/ijod/8157347>
- [6] Li, W., Huang, Z., Cai, R., Yang, W., He, H., & Wang, Y. (2021). Rational Design of Ag/ZnO Hybrid Nanoparticles on Sericin/Agarose Composite Film for Enhanced Antimicrobial Applications. *International Journal of Molecular Sciences*, 22(1), 105. <https://doi.org/10.3390/ijms22010105>